

Synthra Methionine/Choline (Catalog No. 007)

Synthra Methionine/Choline is a completely automated radiosynthesizer for the efficient production of [^{11}C]labeled compounds like [^{11}C]Methionine and [^{11}C]Choline. Automating the synthesis is simple with the easy-to-use configuration software SynthraView. The Synthra Methionine/Choline module offers both, fully automatic as well as manual modes of operation.

Gas Phase Capabilities

- ✓ High specific activities are achieved from in-target produced [^{11}C]CO₂ ranging from 10 Ci/ μmol to 24 Ci/ μmol .*

The target [^{11}C]CO₂ is quantitatively trapped in the stainless steel capillary tubing at -180 °C. After washing out any impurities the CO₂ is released into the methane oven where it is converted to [^{11}C]CH₄ by reduction on a Ni catalyst. Subsequently, the [^{11}C]methane is released and trapped at -120 °C on a Carboxen® CH₄ trap and unreacted hydrogen is removed from the system. In a successive gas phase reaction the [^{11}C]CH₄ is converted into [^{11}C]MeI and trapped on a Porapak Q filled column.


[^{11}C]Labeling Possibilities

- ✓ [^{11}C]Methyl iodide Production: [^{11}C]MeI is ready to be released after 7 minutes starting from the [^{11}C]CO₂ trapping. The yield for the [^{11}C]MeI formation is better than 50% (ndc).
 - Up to **10** sequential methyl iodide preparations are possible from a single box set-up.
 - The [^{11}C]methyl iodide can be used to synthesize [^{11}C]choline by captive chemistry in a loop or for the [^{11}C]Methionine synthesis by solid support heterogeneous reactions on a cartridge.

Additional Gas Phase Options

- ➔ **Methane Option:** A reduced gas phase suitable for the use of CH₄ Target.
- ➔ [^{11}C]HCN (Catalog No.003HCN): The [^{11}C]CH₄ is released with NH₃ gas into a high temperature area where it undergoes a Pt-catalyzed conversion into [^{11}C]HCN at 950 °C.

*Higher specific activities are possible when using methane target.



General Features

- ✓ **Heating and Cooling Capabilities**
 - 6 heating zones
 - 3 with cooling capabilities
 - Temperature range: -196 °C – 1000 °C
- ✓ **Detectors and Controllers**
 - 4 shielded radiation detectors
 - 3 electronic flow controllers (HCN option: 4 flow controller)
 - 1 Pressure sensors as leak detector + 1 filter test unit
- ✓ **Dispensers and Valves**
 - HR-dispenser (up to 50.000 steps, 2.5/5 mL)
 - Chemically inert valves with small dead volume < 35 μL , 5 bar rated

Synthesis Features

- ✓ **4 reagent vials**
 - One small (1 – 3 mL) and three large (10 – 15 mL) volume glass vials for reagents
- ✓ **1 additional Cartridge holders**
- ✓ **SPE Unit** for final product formulation

Synthra C-11 Family

Product Description and Technical Specifications

synthra



✓ Synthesis Options:


- Triflate/Column oven (RT – 200 °C)
- Product Vial Heater Option (Catalog No. 000ph)

Terminal Control

- ✓ **A Laptop (Win 10 Pro) and SynthraView are included**

GMP Features

- ✓ Synthesis files for at least 2 [¹¹C]radiotracers
- ✓ **GMP compliant.** Electronic control and data collection (27/18 channels)
- ✓ **21CFRpart11 & LIMS** compatible



The Graphical User Interface (GUI) of the SynthraView software