

## Synthra HCNplus (Catalog No. 012tv)

Synthra HCNplus is a flexible and completely automated synthesis system for routine production of  $[^{11}\text{C}]\text{HCN}$  on the generation by gas-phase synthesis in combination with a setup for further reaction to the desired tracer product. With easy-to-use configuration software SynthraView, the Synthra HCNplus module offers both, fully automatic and manual modes of operation.

### Gas Phase Capabilities

- ✓ High specific activities are achieved from in-target produced  $[^{11}\text{C}]\text{CO}_2$  ranging from 5 Ci/ $\mu\text{mol}$  to 30 Ci/ $\mu\text{mol}$ .

The target  $[^{11}\text{C}]\text{CO}_2$  is trapped in the stainless steel capillary tubing at  $-180\text{ }^\circ\text{C}$  and converted to  $[^{11}\text{C}]\text{CH}_4$  by reduction on a Ni catalyst. Subsequently, the  $[^{11}\text{C}]\text{methane}$  is trapped at  $-120\text{ }^\circ\text{C}$  on a Carboxen<sup>®</sup>  $\text{CH}_4$  trap and unreacted hydrogen is removed from the system. In a successive gas phase reaction with anhydrous ammonia on a Platinum catalyst at  $950\text{ }^\circ\text{C}$  the  $[^{11}\text{C}]\text{methane}$  is then converted into  $[^{11}\text{C}]\text{HCN}$  and can be used for further synthesis.

### $[^{11}\text{C}]\text{Labeling Possibilities}$

- ✓  **$[^{11}\text{C}]\text{HCN}$  Production:**  $[^{11}\text{C}]\text{HCN}$  is ready for release after only 5 min starting from trapping the  $[^{11}\text{C}]\text{CO}_2$ . The yield is better than 70 %.
  - Up to **50** sequential HCN preparations are possible from a single box set-up.

### Additional Synthesis Options

- ➔ **Methane Option:** A reduced gas phase suitable for the use of  $\text{CH}_4$  Target.
- ➔  **$[^{11}\text{C}]\text{CO}$  (Catalog No. 003CO):** After purification, the  $[^{11}\text{C}]\text{CO}_2$  is released into the column oven for Mo-catalyzed reduction to  $[^{11}\text{C}]\text{CO}$ .

### General Features

- ✓ **Heating and Cooling Capabilities**
  - 6 heating zones
  - 4 with cooling capabilities
  - Temperature range:  $-196\text{ }^\circ\text{C}$  –  $1000\text{ }^\circ\text{C}$



- ✓ **Detectors and Controllers**
  - 6 shielded radiation detectors
  - 4 electronic flow controllers
  - 3 Pressure sensors + 1 filter test unit
- ✓ **Self-Cleaning System**
- ✓ **Dispensers and Valves**
  - HR-dispenser (up to 50.000 steps, 2.5/5 mL)
  - Built-in preparative Radio/UV-HPLC system (0 – 40 mL/min) for product separation and fixed wavelength LED detector with 254 nm or 280 nm
  - 5 spare valves for customization
  - Chemically inert valves with small dead volume < 35  $\mu\text{L}$ , 5 bar rated

### Synthesis Features

- ✓ **2 closed 3 mL reaction vessel** ( $-196\text{ }^\circ\text{C}$  –  $250\text{ }^\circ\text{C}$ ) with integrated cooling to reduce synthesis time (min. volume: 50  $\mu\text{L}$ )
- ✓ **10 reagent vials**
  - Three small (1 – 3 mL) and seven large (10 – 15 mL) volume glass vials for reagents

# Synthra C-11 Family

## Product Description and Technical Specifications

synthra



### ✓ Further Options:


- ➔ Product vial heater option (Catalog No. 000ph)
- ➔ Variable wavelength UV Detector option (Catalog No. 000vuv)

### Terminal Control

- ✓ **A Laptop (Win 10 Pro) and SynthraView are included**

### GMP Features

- ✓ **GMP/GLP compliant.** Electronic control and data collection (27/18 channels)
- ✓ **21CFRpart11 & LIMS** compatible



The Graphical User Interface (GUI) of the SynthraView software